#### Estimado cliente

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Stand der Information: 03/10

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## Dear Customer

Thank you for choosing this high-quality product from Dentaurum. In order to achieve the best results with this product, it is important to carefully study and follow these instructions for use.

The written instructions cannot cover all possible situations that can occur during use. For this reason, our hotline is available to answer any other questions or concerns that may arise.

Due to constant developments, we recommend, even when using the same products regularly, that you study the enclosed current instructions for use or refer to our website at [www.dentaurum.de](http://www.dentaurum.de).

### 1. Manufacturer

Dentaurum GmbH & Co. KG  
Turnstraße 31 · 75228 Ispringen · Germany

### 2. General product description/area of application

The Variety/Variety SP RPE is a universal expansion screw that can be used in the lower arch for creating space, control tipping of the anterior and posterior teeth, as well as opening and closing spaces in both the upper and lower arch. The Variety/Variety SP RPE can be used as a distraction screw, or for early palatal expansion of the upper arch. The screw is welded or soldered to preformed bands using its retention legs and does not require the use of acrylic. The bands are then cemented into place.

### Variety (Fig. 1)

- small dimensions with maximum expansion
- solder free
- break resistant retention legs and stable guidance throughout the entire expansion
- the integrated spindle friction prevents unwanted turning back
- permanently laser marked arrow on the body of the screw indicates the opening direction and a figure the maximum expansion
- screw can be tied down once desired expansion has been reached

### Variety SP (small palate, Fig. 2)

- all the characteristics of the Variety expansion screw
- especially suitable for very narrow jaws
- makes it possible to bend the retention legs very close to the screw body

### 3. Fabricating a Variety/Variety SP appliance

Preformed bands are fitted on the mandibular left and right 1st bicuspid and left and right 1st molars on a hard plaster model. The legs of the Variety/Variety SP screw are then bent to fit the bands so that they can be soldered or welded into place (Fig. 3).

The Variety/Variety SP screw should „float“ freely 1-2 mm above the gingiva (Fig. 4). To avoid any leverage effect, please be sure the direction of force runs parallel to the occlusal plane.

### 4. Working tips for the Variety/Variety SP screws

In order to bend the retention legs of the Variety/Variety SP screw gently and easily, we recommend using the universal bending key (REF 611-122-00, Fig. 5 + 9) without the aid of pliers (lever principle). If this is not available, the retention legs must be stabilised by holding at the expansion screw body with a pair of sturdy pliers e.g. Wire bending pliers Maxi (REF 013-522-00), in order to prevent the retention legs from breaking (Fig. 6).

When using the soldering technique please protect the body of the screw against excessive heat, e.g. with heat screen (REF 380-070-00).

### 5. Note for doctor and patient

Before the Variety/Variety SP appliance is cemented into position, turn the spindle (screw) slightly (2-3 quarter turns) in the direction of the arrow to test the mobility. Then turn back to the original position (Fig. 7).

The finished Variety/Variety SP appliance is cemented in place in the conventional manner. The cement manufacturer's instructions must be observed.

Following this, the screw is turned four times by 90° using the key (REF 611-120-01, Fig. 8 + 10) supplied for this purpose. One full 360° turn of the spindle = 0.8 mm expansion.

The patient should turn the screw two times 90° both in the morning and evening (according to your instructions) until the desired turn is achieved.

**Important:** use only the blue safety key (Fig. 10) supplied, as this cannot accidentally be swallowed.

The safety key is designed to hold a maximum force of 35 N  $\pm$  3.5 kg. Improper use may lead to overload and consequently to failure thus potentially excluding any guarantee and liability.

Regular visits are recommended to monitor progress.

The Variety/Variety SP appliance is secured against unwanted turning-back, after the desired degree of expansion is achieved, by fitting a ligature or blocking with acrylic. For stabilization the appliance should remain in the mouth for approximately 2-3 months.

A retainer should be worn following removal of Variety/Variety SP appliance.

### 6. Composition

The composition is included in the materials list, refer to the catalog or [www.dentaurum.de](http://www.dentaurum.de).

### 7. Safety information

Always use the blue safety key supplied to prevent swallowing!

The product should not be used, if there is a known allergic reaction to one or more of the material components.

### 8. Information on single-use products

The product is intended for single use only. Reconditioning of a product that has already been used (recycling) and its reuse on a patient is not permitted.

### 9. Availability

Description	REF	Quantity
Variety	REF 602-806-10	1 piece
Variety	REF 602-806-30	10 pieces
Variety SP	REF 602-809-10	1 piece
Variety SP	REF 602-809-30	10 pieces
with safety key		
Safety key	REF 611-120-01	1 piece
Universal bending key	REF 611-122-00	1 piece

### 10. With reference to our quality standards

Dentaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

### 11. Explanation of symbols used on the label

REF	Order number
	Batch description
0483	CE Mark MP Class IIa
<a href="http://www.dentaurum.de">www.dentaurum.de</a>	Instructions for use available in the internet
	Not for reuse
	Codablock barcode
	Matrix code
Date of information: 03/10	
Subject to modifications	

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